



# Applied Scientific Instrumentation

## Stage Inserts

### 160 x 110 mm PETRI DISH and FLASK Inserts

#### Universal Petri Dish Insert

I-3026

The I-3026 insert accepts 30mm to 65 mm (1.18" to 2.56") Petri dishes. The unit places the bottom of a dish about 7.4 mm below the top of the insert.

Depth from Top of Insert: 7.4 mm  
Overall Thickness: 8.0 mm



#### 39 mm Petri Dish Insert

I-3039

The I-3039 insert accepts Petri dishes with a maximum bottom outside diameter of 39.0 mm (1.54"). The unit places the bottom of the dish about 7.4 mm below the top of the insert. The surrounding cavity diameter is 46.6 mm (1.83") to accommodate a dish cover. Hole diameter is 36.4 mm (1.43"). Ease of access side slots. Spring clips for secure stability.

Depth from Top of Insert: 7.4 mm  
Overall Thickness: 8.0 mm



#### 35 mm & 50 mm WillCo-dish Insert

I-3034

The I-3034 insert accepts WillCo dishes with maximum bottom outside diameters of 33.9 mm (1.33") and 50.0 mm (1.97"). The unit places the bottom of the dishes about 7.3 mm below the top of the insert. Surrounding cavity diameters are 37.7 mm (1.48") and 54.0 mm (2.12") to accommodate dish covers. Hole diameters are 31.3 mm (1.23") and 47.5 mm (1.97"). Ease of access side slots. Spring clips for secure stability.

Depth from Top of Insert: 7.3 mm  
Overall Thickness: 8.0 mm



#### 50 mm Petri Dish Insert

I-3050

The I-3050 insert accepts Petri dishes with a maximum bottom outside diameter of 50.0 mm (1.97"). The unit places the bottom of the dish about 7.4 mm below the top of the insert. The surrounding cavity diameter is 53.9 mm (2.12") to accommodate a dish cover. Hole diameter is 47.5 mm (1.87"). Ease of access side slots. Spring clips for secure stability.

Depth from Top of Insert: 7.4 mm  
Overall Thickness: 8.0 mm



#### 35 mm Petri Dish Insert

I-3035

The I-3035 insert accepts Petri dishes with a maximum bottom outside diameter of 35.3 mm (1.39"). The unit places the bottom of the dish about 7.4 mm below the top of the insert. The surrounding cavity diameter is 41.2 mm (1.62") to accommodate a dish cover. Hole diameter is 31.8 mm (1.25"). Ease of access side slots. Spring clips for secure stability.

Depth from Top of Insert: 7.4 mm  
Overall Thickness: 8.0 mm




#### 60 mm Petri Dish Insert

I-3060

The I-3060 insert accepts Petri dishes with a maximum outside diameter of 61.9 mm (2.44"). The unit places the bottom of the dish about 7.4 mm below the top of the insert. Surrounding cavity diameter is 67.0 mm (2.64") to accommodate a dish cover. Hole diameter is 57.0 mm (2.25"). Ease of access side slots. Spring clips for secure stability.

Depth from Top of Insert: 7.4 mm  
Overall Thickness: 8.0 mm



*We Create Solutions*

Applied Scientific Instrumentation, Inc. ♦ 29391 W Enid Rd ♦ Eugene, OR USA 97402-9533  
(541) 461-8181 ♦ (800) 706-2284 ♦ info@ASImaging.com ♦ www.ASImaging.com

### 85 mm Petri Dish Insert

**I-3085**

The I-3085 insert accepts Petri dishes with a maximum bottom outside diameter of 88.9 mm (3.50"). The unit places the bottom of the dish about 7.4 mm below the top of the insert. Hole diameter is 82.5 mm (3.25"). Ease of access side slots. Spring clips for secure stability.

Depth from Top of Insert: 7.4 mm  
Overall Thickness: 8.0 mm



### 88 Petri Dish Insert

**I-3088**

The I-3088 insert accepts Petri dishes with a maximum bottom outside diameter of 90.9 mm (3.58"). The unit places the bottom of the dish about 7.4 mm below the top of the insert. Hole diameter is 84.2 mm (3.31"). Ease of access side slots. Spring clips for secure stability.

Depth from Top of Insert: 7.4 mm  
Overall Thickness: 8.0 mm



### 98 mm Petri Dish Insert

**I-3098**

The I-3098 insert accepts Petri dishes with a maximum bottom outside diameter of 100.8 mm (3.97"). The unit places the bottom of the dish about 7.4 mm below the top of the insert. Hole diameter is 94.5 mm (3.72"). Ease of access side slots. Spring clips for secure stability.

Depth from Top of Insert: 7.4 mm  
Overall Thickness: 8.0 mm



### Flask Insert for Corning T25 Flask

**I-3090-C**

The I-3090-C insert accepts Corning T25 flasks. The unit places the bottom of the flask about 7.3 mm below the top of the insert, and uses a spring-loaded clamp for stability.

Depth from Top of Insert: 7.3 mm  
Overall Thickness: 8.0 mm



### Flask Insert for Corning T75 Flask

**I-3092-C**

The I-3092-C insert accepts Corning T75 flasks. The unit places the bottom of the flask about 7.3 mm below the top of the insert, and uses a spring-loaded clamp for stability. The I-3092-C also holds SARSTEDT-brand T75 flasks.

Depth from Top of Insert: 7.3 mm  
Overall Thickness: 8.0 mm




### Flask Insert for Nunc T25 Flask

**I-3090-N**

The I-3090-N insert accepts Nunc T25 flasks. The unit places the bottom of the flask about 7.3 mm below the top of the insert, and uses a spring-loaded clamp for stability.

Depth from Top of Insert: 7.3 mm  
Overall Thickness: 8.0 mm



### Flask Insert for Nunc T75 Flask

**I-3092-N**

The I-3092-N insert accepts Nunc T75 flasks. The unit places the bottom of the flask about 7.3 mm below the top of the insert, and uses a spring-loaded clamp for stability.

Depth from Top of Insert: 7.3 mm  
Overall Thickness: 8.0 mm



***We Create Solutions***

Applied Scientific Instrumentation, Inc. ♦ 29391 W Enid Rd ♦ Eugene, OR USA 97402-9533  
(541) 461-8181 ♦ (800) 706-2284 ♦ info@ASImaging.com ♦ www.ASImaging.com